



TAG-72


USA IV D

This is directed against a high molecular weight glycoprotein known as tumor-associated glycoprotein-72.

Mouse Monoclonal, Clone: B72.3

Control: Lung Adenocarcinoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

Shelf life: Up to 36 months. **Isotype:** IgG1/K **Cellular localization:** Cytoplasmic



Lung Adenocarcinoma

Catalog no.

TA-2010-1

TA-2010-2

TA-2010-3

TA-2010-4

TA-2010-5

TA-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides