




# TRAcP

USA IVD

This labels the cells of hairy cell leukemia with high sensitivity and specificity.

**Mouse Monoclonal, Clone: 9C5**

**Control:** Hairy Cell Leukemia **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>2b</sub> **Cellular localization:** Cytoplasmic



Hairy Cell Leukemia

## Catalog no.

TR-2010-1

TR-2010-2

TR-2010-3

TR-2010-4

TR-2010-5

TR-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides