




# TSH

USA IV D

This is a marker in classification of pituitary tumors and differential identification of primary and metastatic tumors in the pituitary gland.

## Rabbit Polyclonal

**Control:** Pituitary **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Pituitary

## Catalog no.

TS-2010-1

TS-2010-2

TS-2010-3

TS-2010-4

TS-2010-5

TS-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides