



# TTF-1


USA IVD

This is useful in differentiating primary adenocarcinoma of the lung from metastatic carcinomas from the breast, mediastinal germ cell tumors and malignant mesothelioma.

**Mouse Monoclonal, Clone: 8G7G3/1**

**Control:** Lung Adenocarcinoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

**Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Nuclear



Lung Adenocarcinoma

## Catalog no.

TT-2010-1

TT-2010-2

TT-2010-3

TT-2010-4

TT-2010-5

TT-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides