




TdT

USA IVD

This is useful for the diagnosis of T-ALL/T-lymphoblastic lymphoma and B-ALL/B-lymphoblastic lymphoma.

Rabbit Polyclonal

Control: Thymus **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Cellular localization:** Nuclear



Thymus

Catalog no.

TD-2010-1

TD-2010-2

TD-2010-3

TD-2010-4

TD-2010-5

TD-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides