




Thyroglobulin

USA IVD

This can identify medullary carcinomas of the thyroid and can differentiate between primary thyroid and lung neoplasms.

Mouse Monoclonal, Clone: 2H11+6E1

Control: Thyroid **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic



Thyroid

Catalog no.

TG-2010-1

TG-2010-2

TG-2010-3

TG-2010-4

TG-2010-5

TG-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides