




Tyrosinase

USA IVD

This is specific for malignant melanoma and melanotic neurofibroma.

Mouse Monoclonal, Clone: T311

Control: Melanoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2a} **Cellular localization:** Cytoplasmic



Melanoma

Catalog no.

TY-210-1

TY-2010-2

TY-2010-3

TY-2010-4

TY-2010-5

TY-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides