




Villin

USA IVD

This is useful in differentiating gastrointestinal adenocarcinomas, neuroendocrine carcinomas, and ovarian adenocarcinomas from others.

Mouse Monoclonal, Clone: CWWB1

Control: Colon **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic, Membranous



Colon

Catalog no.

VI-2010-1

VI-2010-2

VI-2010-3

VI-2010-4

VI-2010-5

VI-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.