




Vimentin

USA IVD

This can distinguish melanomas from undifferentiated carcinomas and large cell lymphomas.

Mouse Monoclonal, Clone: V9

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG/K **Cellular localization:** Cytoplasmic



Tonsil

Catalog no.

VM-2010-1

VM-2010-2

VM-2010-3

VM-2010-4

VM-2010-5

VM-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides