

# CYTOLOGY AND HISTOLOGY SUPPLIES

## PROCESSING / EMBEDDING TISSUE CASSETTES



Same capacity and displacement as other common tissue cassettes.

Covers hold specimens in place but may be opened and closed as often as needed. Covers may be detached and mated to new coverless cassettes to reduce both costs and plastic waste. Available in white, blue, green, gray, pink, yellow, orange, light brown and red.

### CATALOG #

### Description

H-4908      Processing / embedding tissue cassette with covers



Multi-compartment Biopsy Cassettes w/ covers (2, 4 and 9 compartments). Colors: white, green and gray.

### CATALOG #

### Description

H-12408-2      Multi-compartments cassette w/ covers

H-12408-4

H-12408-9

**Available as packs of 50, 100, 250 or 500**

# CYTOLOGY AND HISTOLOGY SUPPLIES

## COPLIN STAINING JARS AND DISHES



Coplin staining jars. The cover and the body are made of heavy-duty molded glass. Holds 5 slides. Slots are wide enough to hold 10 slides (back-to-back) Height 4-1/4" (108 mm) and the O.D. of the cover is 2-3/8" (60mm)

CATALOG #

H-221108

Description

Coplin staining jars



Dimensions of 3.75"x3" x1.75" (95 mm x 76 mm x 44 mm)

CATALOG #

H-231108

Description

Coplin staining dishes