

CYTOLOGY AND HISTOLOGY SUPPLIES

FORCEPS



Self-sprung. Transparent with blue tint. (Sterile and non-sterile)



Self-sprung. Green opaque. Constructed of acetal. (Resistant to acid and alkali)



Self-sprung with pointed end. Available in blue only. (Sterile and non-sterile)

<u>CATALOG #</u>	<u>Description</u>	<u>Length</u>
H-5908	Blue plastic forceps	130mm
H-5908S	Blue plastic forceps (Individual sterile)	130mm
H-12208	Green plastic forceps	105mm
H-12308	Transparent with blue tint	105mm
H-12308S	Transparent with blue tint (Individual sterile)	105mm

CYTOLOGY AND HISTOLOGY SUPPLIES

SLIDE STORAGE / MAILING BOXES



Heavy duty slide storage boxes. New slide storage boxes are made of plastic for long lasting wear. They have user-friendly hinge with an easy to use latch. They are stackable for space saving. Available in blue, cream, green, yellow, white, and red.

CATALOG

Description

H-1908

25 slide box

H-2908

50 slide box

H-3908

100 slide box



SLIDE STORAGE / MAILING BOXES

Made of a styrenic thermoplastic material.

CATALOG

Description

H-12508-S/M

25 slide box

Sold as packs of 100