




# p21 WAF1

USA IVD

This is a nuclear 21kD protein and a product of WAF1/CIP1 gene.

**Mouse Monoclonal, Clone: DCS-60.2**

**Control:** Colon **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>2a</sub> **Cellular localization:** Nuclear



Colon

## Catalog no.

PW-2010-1

PW-2010-2

PW-2010-3

PW-2010-4

PW-2010-5

PW-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides