




p27 Kip1

USA IVD

This is a human gene which encodes a protein which belongs to the Cip/Kip family of cyclin dependent kinase inhibitor proteins.

Mouse Monoclonal, Clone: SX53G8

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Nuclear



Tonsil

Catalog no.

PK-2010-1

PK-2010-2

PK-2010-3

PK-2010-4

PK-2010-5

PK-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .