




## p57 Kip2

USA IVD

This is used in the identification of complete hydatidiform mole from partial hydatidiform mole and hydropic abortion.

**Mouse Monoclonal, Clone: Kp10**

**Control:** Placenta **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG2b/K **Cellular localization:** Nuclear



Placenta

### Catalog no.

PI-2010-1

PI-2010-2

PI-2010-3

PI-2010-4

PI-2010-5

PI-2010-6

sizes available upon  
request

### Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

3 ml Prediluted

7 ml Prediluted

5 Positive control slides

**Prediluted.** Smallest