




## p63

USA RUO

This is used to differentiate benign versus malignant prostate and breast lesions.

**Mouse Monoclonal, Clone: 4A4**

**Control:** Prostate **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8<sup>0</sup>C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>2a</sub>/K **Cellular localization:** Nuclear



Prostate

### Catalog no.

PS-2010-1

PS-2010-2

PS-2010-3

PS-2010-4

PS-2010-5

PS-2010

### Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides